

Product Datasheet

ADH7 antibody (orb624518)

Description	Rabbit polyclonal antibody to ADH7
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	alcohol dehydrogenase 7(class IV), mu or sigma polypeptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	41 kDa, observed 45-47 kDa
Uniprot ID	P40394
Dilution Range	WB: 1:500-1:5000, IF: 1:20-1:200, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.