

Product Datasheet

ADIPOQ antibody (orb624522)

Description	Mouse monoclonal antibody to ADIPOQ
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	adiponectin, C1Q and collagen domain containing
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	26 kDa
Uniprot ID	Q15848
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.