

Product Datasheet

AGPAT9 antibody (orb624557)

Description	Rabbit polyclonal antibody to AGPAT9
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	1-acylglycerol-3-phosphate O-acyltransferase 9
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	49 kDa
Uniprot ID	Q53EU6
Dilution Range	WB: 1:500-1:5000, IF: 1:20-1:200, IHC-P: 1:50- 1:400
Expiration Date	12 months from date of receipt.