

Product Datasheet

AKAP14 antibody (orb624587)

Description	Rabbit polyclonal antibody to AKAP14
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	A kinase(PRKA) anchor protein 14
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	23 kDa, observed 28 kDa
Uniprot ID	Q86UN6
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:20-200
Expiration Date	12 months from date of receipt.