

Product Datasheet

ALKBH3 antibody (orb624641)

Description	Rabbit polyclonal antibody to ALKBH3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	alkB, alkylation repair homolog 3(E. coli)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	33 kDa, observed 25 kDa
Uniprot ID	Q96Q83
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.