

Product Datasheet

Alpha actin antibody (orb624647)

Description	Rabbit polyclonal antibody to Alpha actin
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	actin, alpha, cardiac muscle 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	42 kDa
Uniprot ID	P68032
Dilution Range	WB: 1:200-1:2000, IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.