

Product Datasheet

ANKS1B antibody (orb624725)

Description	Rabbit polyclonal antibody to ANKS1B
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	ankyrin repeat and sterile alpha motif domain containing 1B
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	138 kDa, Observed 70 kDa
Uniprot ID	Q7Z6G8
Dilution Range	WB: 1:500-1:5000, IP: 1:500-1:5000, IHC-P: 1:20-1:200, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.