

Product Datasheet

APEX1 antibody (orb624772)

Description	Rabbit polyclonal antibody to APEX1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	APEX nuclease(multifunctional DNA repair enzyme) 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	36 kDa
Uniprot ID	P27695
Dilution Range	WB: 1:500 - 1:2000; IHC: 1:50 - 1:100; IF: 1:50 - 1:100
Expiration Date	12 months from date of receipt.