

Product Datasheet

APOBEC3G antibody (orb624787)

Description	Rabbit polyclonal antibody to APOBEC3G
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	apolipoprotein B mRNA editing enzyme, catalytic polypeptide-like 3G
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	46 kDa
Uniprot ID	Q9HC16
Dilution Range	WB: 1:200-1:1000; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.