

Product Datasheet

ARMC5 antibody (orb624864)

Description	Rabbit polyclonal antibody to ARMC5
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	armadillo repeat containing 5
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	55 kDa
Uniprot ID	Q96C12
Dilution Range	WB: 1:200-1:1000
Expiration Date	12 months from date of receipt.