

Product Datasheet

ATG2B antibody (orb624941)

Description	Rabbit polyclonal antibody to ATG2B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	ATG2 autophagy related 2 homolog B(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	240-250 kDa
Uniprot ID	Q96BY7
Dilution Range	WB: 1:500-1:1000, IHC-P: 1:50-1:500, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.