



Product Datasheet

ITLN1 antibody (orb6250)


Description	ITLN1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ITLN1 (131-230/313 aa)
Target	ITLN1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35 kDa
Uniprot ID	Q8WWAO
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse embryo...



WB analysis of mouse heart Lysate at 40 ...



WB analysis of human MCF-7 cell Lysate a...