

Product Datasheet

B-Myb antibody (orb625157)

Description	Rabbit polyclonal antibody to B-Myb
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	v-myb myeloblastosis viral oncogene homolog(avian)-like 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	90 kDa
Uniprot ID	P10244
Dilution Range	WB: 1:500-1:2000; IP: 1:500-1:1000; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.