

Product Datasheet

BPNT1 antibody (orb625168)

Description	Rabbit polyclonal antibody to BPNT1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	3'(2'), 5'-bisphosphate nucleotidase 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	36-40 kDa
Uniprot ID	O95861
Dilution Range	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:500
Expiration Date	12 months from date of receipt.