



Product Datasheet

Jak3 antibody (orb6252)


Description	Jak3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Jak3 (601-700/1124 aa)
Target	Jak3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125 kDa
Uniprot ID	P52333
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=1ug/test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Flow cytometry of Molt-4.Primary Antibod...



Western blot:SP2/0 Cell(Mouse);Jurkat Ce...



Western blot:Lane 1:Mouse Thymus;Lane 2:...