

Product Datasheet

C11orf57 antibody (orb625223)

Description	Rabbit polyclonal antibody to C11orf57
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	chromosome 11 open reading frame 57
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30-34 kDa
Uniprot ID	Q6ZUT1
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.