

Product Datasheet

C19orf61 antibody (orb625252)

Description	Rabbit polyclonal antibody to C19orf61
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	chromosome 19 open reading frame 61
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	57-65 kDa
Uniprot ID	Q9H0W8
Dilution Range	WB: 1:500-1:2000; IP: 1:500-1:1000; IHC: 1:20-1:200; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.