

Product Datasheet

C1R antibody (orb625268)

Description	Rabbit polyclonal antibody to C1R
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	complement component 1, r subcomponent
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	85 kDa, 30 kDa
Uniprot ID	P00736
Dilution Range	WB: 1:500-1:2000; IF: 1:50-1:200
Expiration Date	12 months from date of receipt.