

Product Datasheet

Caspase 9 antibody (orb625441)

Description	Rabbit polyclonal antibody to Caspase 9
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	caspase 9, apoptosis-related cysteine peptidase
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	46kd
Uniprot ID	P55211
Dilution Range	WB: 1:200-1:1000, IP: 1:200-1:1000, IHC-P: 1:50-1:200, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.