

Product Datasheet

CC2D1B antibody (orb625479)

Description	Rabbit polyclonal antibody to CC2D1B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	coiled-coil and C2 domain containing 1B
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	120 kDa
Uniprot ID	Q5T0F9
Dilution Range	WB: 1:500-1:1000; IP: 1:200-1:1000
Expiration Date	12 months from date of receipt.