

## Product Datasheet

### CD109 antibody (orb625542)

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to CD109                                                                                   |
| <b>Species/Host</b>        | Rabbit                                                                                                                |
| <b>Reactivity</b>          | Human                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Tested Applications</b> | ELISA, WB                                                                                                             |
| <b>Immunogen</b>           | CD109 molecule                                                                                                        |
| <b>Form/Appearance</b>     | liquid                                                                                                                |
| <b>Storage</b>             | PBS with 0.02% sodium azide and 50% glycerol<br>pH 7.3, -20°C for 24 months (Avoid repeated<br>freeze / thaw cycles.) |
| <b>Note</b>                | For research use only                                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                            |
| <b>Purity</b>              | ≥95% as determined by SDS-PAGE                                                                                        |
| <b>MW</b>                  | 180 kDa, 150 kDa, 120 kDa, 75 kDa                                                                                     |
| <b>Uniprot ID</b>          | <b>Q6YHK3</b>                                                                                                         |
| <b>Dilution Range</b>      | WB: 1:500-1:2000                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                       |