

Product Datasheet

CD11a/Integrin alpha-L antibody (orb625543)

Description	Mouse monoclonal antibody to CD11a/Integrin alpha-L
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	integrin, alpha L(antigen CD11A(p180), lymphocyte function-associated antigen 1; alpha polypeptide)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	129 kDa
Uniprot ID	P20701
Dilution Range	IHC: 1:20-1:500
Expiration Date	12 months from date of receipt.