

Product Datasheet

CD138/Syndecan-1 antibody (orb625556)

Description	Mouse monoclonal antibody to CD138/Syndecan-1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	syndecan 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	60-100 kDa
Uniprot ID	P18827
Dilution Range	WB: 1:500-1:1000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.