

Product Datasheet

Chromogranin A antibody (orb625796)

Description	Rabbit polyclonal antibody to Chromogranin A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	chromogranin A(parathyroid secretory protein 1)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa, 75 kDa
Uniprot ID	P10645
Dilution Range	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.