

Product Datasheet

COX10 antibody (orb625963)

Description	Rabbit polyclonal antibody to COX10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	COX10 homolog, cytochrome c oxidase assembly protein, heme A: farnesyltransferase(yeast)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	51 kDa
Uniprot ID	Q12887
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.