

Product Datasheet

COX IV antibody (orb625985)

Description	Mouse monoclonal antibody to COX IV
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, WB
Immunogen	cytochrome c oxidase subunit IV isoform 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	17-18 kDa
Uniprot ID	P13073
Dilution Range	WB: 1:2000-1:4000; IHC: 1:200-1:500; IF: 1:50- 1:500
Expiration Date	12 months from date of receipt.