

Product Datasheet

CREB5 antibody (orb626025)

Description	Rabbit polyclonal antibody to CREB5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	cAMP responsive element binding protein 5
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60-70 kDa
Uniprot ID	Q02930
Dilution Range	WB: 1:500-1:3000
Expiration Date	12 months from date of receipt.