

Product Datasheet

CSTF3 antibody (orb626090)

Description	Rabbit polyclonal antibody to CSTF3
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	cleavage stimulation factor, 3' pre-RNA, subunit 3, 77kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	77-83 kDa
Uniprot ID	Q12996
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.