

Product Datasheet

CWC25 antibody (orb626132)

Description	Rabbit polyclonal antibody to CWC25
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	coiled-coil domain containing 49
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	Q9NxE8
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.