

Product Datasheet

CXCL17 antibody (orb626142)

Description	Rabbit polyclonal antibody to CXCL17
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	chemokine(C-X-C motif) ligand 17
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q6UXB2
Dilution Range	IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.