

Product Datasheet

Cytokeratin 20 antibody (orb626235)

Description	Mouse monoclonal antibody to Cytokeratin 20
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	keratin 20
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	46 kDa
Uniprot ID	P35900
Dilution Range	WB: 1:500-1:2000; IHC: 1:1000-1:10000
Expiration Date	12 months from date of receipt.