

Product Datasheet

DCLRE1B antibody (orb626280)

Description	Rabbit polyclonal antibody to DCLRE1B
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	DNA cross-link repair 1B(PSO2 homolog, S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa
Uniprot ID	Q9H816
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.