

Product Datasheet

DDB1 antibody (orb626298)

Description	Mouse monoclonal antibody to DDB1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	damage-specific DNA binding protein 1, 127kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	127 kDa
Uniprot ID	Q16531
Dilution Range	WB: 1:500-1:2000; IHC: 1:200-1:800
Expiration Date	12 months from date of receipt.