

Product Datasheet

DNAJB1 antibody (orb626437)

Description	Rabbit polyclonal antibody to DNAJB1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	DnaJ(Hsp40) homolog, subfamily B, member 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	40kd
Uniprot ID	P25685
Dilution Range	WB: 1:500-1:3200, IP: 1:500-1:1000, IHC-P: 1:20-1:200, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.