

Product Datasheet

DNAJC15 antibody (orb626453)

Description	Rabbit polyclonal antibody to DNAJC15
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	DnaJ(Hsp40) homolog, subfamily C, member 15
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	18-20kd
Uniprot ID	Q9Y5T4
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.