

Product Datasheet

DOCK3 antibody (orb626485)

Description	Rabbit polyclonal antibody to DOCK3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	dedicator of cytokinesis 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	233 kDa
Uniprot ID	Q8IZD9
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.