

Product Datasheet

Estrogen Receptor beta antibody (orb626823)

Description	Mouse monoclonal antibody to Estrogen Receptor beta
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	estrogen receptor 2(ER beta)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	50 kDa
Uniprot ID	Q92731
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.