


Product Datasheet

KIF17 antibody (orb6269)

Description	KIF17 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human KIF17 (201-300/1029 aa)
Target	KIF17
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	115 kDa
Uniprot ID	Q9P2E2
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 lysates u...