

Product Datasheet

Flag tag antibody (orb627073)

Description	Rabbit polyclonal antibody to Flag tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Flag tag
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:1000-1:4000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.