

Product Datasheet

Galectin 8 antibody (orb627209)

Description	Rabbit polyclonal antibody to Galectin 8
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	lectin, galactoside-binding, soluble, 8
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	000214
Expiration Date	12 months from date of receipt.