

Product Datasheet

GAR1 antibody (orb627238)

Description	Rabbit polyclonal antibody to GAR1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	GAR1 ribonucleoprotein homolog(yeast)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	28kd
Uniprot ID	Q9NY12
Dilution Range	WB: 1:500-1:1000; IP: 1:200-1:1000; IHC: 1:20- 1:200; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.