

Product Datasheet

Glutamate receptor 2 antibody (orb627379)

Description	Rabbit polyclonal antibody to Glutamate receptor 2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	glutamate receptor, ionotropic, AMPA 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	94-105 kDa
Uniprot ID	P42262
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.