

Product Datasheet

GPR50 antibody (orb627476)

Description	Rabbit polyclonal antibody to GPR50
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC
Immunogen	G protein-coupled receptor 50
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	67 kDa
Uniprot ID	Q13585
Dilution Range	IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.