

Product Datasheet

GSTO1 antibody (orb627551)

Description	Mouse monoclonal antibody to GSTO1
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	glutathione S-transferase omega 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	28kd
Uniprot ID	P78417
Dilution Range	WB: 1:1000-1:4000; IHC: 1:50-1:200
Expiration Date	12 months from date of receipt.