

Product Datasheet

HMGCS1 antibody (orb627742)

Description	Rabbit polyclonal antibody to HMGCS1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	3-hydroxy-3-methylglutaryl-Coenzyme A synthase 1(soluble)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	57 kDa
Uniprot ID	Q01581
Dilution Range	WB: 1:1000-1:10000, IP: 1:500-1:5000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.