

Product Datasheet

HPRT1 antibody (orb627801)

Description	Rabbit polyclonal antibody to HPRT1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	hypoxanthine phosphoribosyltransferase 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	24-28 kDa
Uniprot ID	P00492
Dilution Range	WB: 1:1000-1:10000, IP: 1:500-1:5000, IHC-P: 1:20-1:200, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.