

Product Datasheet

IA-2 antibody (orb627886)

Description	Mouse monoclonal antibody to IA-2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	protein tyrosine phosphatase, receptor type, N
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	60-65kd, 106-115kd
Uniprot ID	Q16849
Dilution Range	WB: 1:500-1:1000; IHC: 1:200-1:800; IF: 1:10-1:100
Expiration Date	12 months from date of receipt.