



Product Datasheet

Interferon gamma antibody (orb627889)


Description	Rabbit polyclonal antibody to Interferon gamma
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	allograft inflammatory factor 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	17 kDa
Uniprot ID	P55008
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:50-1:500
Expiration Date	12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Immunofluorescent analysis of (-20°C E...



Various lysates were subjected to SDS PA...