

Product Datasheet

IFIT3 antibody (orb627924)

Description	Rabbit polyclonal antibody to IFIT3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	interferon-induced protein with tetratricopeptide repeats 3
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60 kDa
Uniprot ID	O14879
Dilution Range	WB: 1:500-1:1000; IP: 1:500-1:1000
Expiration Date	12 months from date of receipt.