

Product Datasheet

IFT27 antibody (orb627938)

Description	Rabbit polyclonal antibody to IFT27
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	RAB, member of RAS oncogene family-like 4
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	20 kDa
Uniprot ID	Q9BW83
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.